













SEQUIRE®

Microcatheter / Mikrokatetr / Mikrokateter / Μικροκαθετήρας / Microcateter / Mikrokatéter / Microcathéter / Microcatetere / Microcatéter / Mikrocewnik / Бичил катетер / สายสวนหลอดเลือดเล็กพิเศษ / Ống thông siêu nhỏ / 微導管 / 마이크로카테터

ENGLISH Instructions for use	5	NORSK Bruksanvisning	47
Français Notice d'utilisation	8	ROMÂNĂ Instrucțiuni de utilizare	50
Deutsch Gebrauchsanweisung	12	SLOVENSKY Návod na použitie	53
Italiano Istruzioni per l'uso	15	SVENSKA Användningsinstruktioner	57
Español Instrucciones de uso	19	Magyar Használati utasítás	60
Nederlands Gebruiksaanwijsties	22	Bahasa Melayu Arahan penggunaan	64
Português Instruções de utilização	26	Монгол Хэл Хэрэглэх Заавар	67
Polski Instrukcja Obsługi	30	ไทย คำแนะนำในการใช้งาน	71
Türkçe Kullanım Kılavuzu	33	Tiếng Việt Hướng dẫn sử dụng	74
ČESKY Návod k použití	36	漢語 使用說明	77
DANSK Brugervejledning	40	한국어 사용 설명서	80
ΕΛΛΗΝΙΚΑ Οδηγίες χρήσης	43		

Compatible embolic bead (microsphere) size/ Taille de Billes (microsphères) d'embolisation compatible/ Kompatible Mikrosphärengröße (Embolisation)/ Formato microsfere embolizante compatibile/

Tamaño compatible de la esfera (microesfera) embolizante/ Compatibele grootte embolisch materiaal (microfeer)/ Tamanho de microesferas de embolização compatíveis/ Zgodna wielkość kulek (mikrosfer) do embolizacji/ Uygun emboilk küre (mikrosfer) boyutu/ Velikost kompatibilnih emboličkih koralū (mikrosfery)/ Kompatibel størrelse af emboliske partikler (mikrosfærer)/ Μέγεθος συμβατού σφαιριδίου (μικροσφαιριδίου) εμβολισμού/ Kompatibel mikrosfærørørrelse (embolisering)/ Dimensiune compatibilă cu particulele de embolizare (microsfere)/ Kompatibilná veľkosť embolizačných guľičiek (mikrosfery)/ Kompatibel embolsk Bead Storlek/ Kompatibilis embolizációs gyöngy (mikrogömb) mérete/ Saiz bebola (mikrosfera) embolisk yang sesuai/ Тохирхэм эмбол дүслүн (бичил дугиар) хэмжээ/ ขนาดเม็ดฝัง (ไมโครสเฟียร์) ที่ใช้ได้/ Kích thước hạt (vi cầu) gây thromb tắc tương thích/ 相容的栓塞球 (微球體) 大小/ 호환되는 색전 구슬 (미크로구체) 크기

Please consult the accompanying documentation/ Prière de consulter la documentation jointe/ Begleitedokumente lesen/ Consultare la documentazione di accompagnamento/ Consulte la documentación adjunta/ Raadpleeg de begeleidende documentatie/ Consulte os documentos anexos/ Zapoznać się z dokumentacją towarzyszącą/ Lüften ekte verilen belgelere bakın/ Prosim pročítate si doprovodné dokumenty/ Se veligst den medfølgende dokumentation/ Συμβουλευτείτε τα συνοδευτικά έγγραφα/ Se medfølgende dokumentasjon/ Vá rugám szá konzultálni dokumentációja anexát/ Prečítajte si sprievodnú dokumentáciu/ Se bifogad dokumentation/ Olvassa el a kísérődokumentációt/ Sila rukuj kepada dokumen yang disertakan dengan peranti ini/ Хабарсан бичиг баримтай танилцана уу/ โปรดศึกษาเอกสารที่แนบมาด้วย/ Vui lòng tham khảo tài liệu kèm theo/ 請參考隨附文件/ 함께 제공되는 설명서를 참고하십시오

Caution! Consult accompanying documentation/ Attention ! Consulter la documentation jointe/ Achtung! Begleitedokumente lesen/ Attenzione! Consultare la documentazione di accompagnamento/ ¡Precaución! Consulte la documentación adjunta/ Let op! Raadpleeg de begeleidende documentatie/ Recomendaciones! Consulte os documentos anexos/ Przestoga! Zapoznać się z dokumentacją towarzyszącą/ Dikkat! Ekte verilen belgelere bakın/ Pozor! Pročítate si doprovodné dokumenty/ Forsigtig! Se den medfølgende dokumentation/ Προσοχή! Συμβουλευτείτε τα συνοδευτικά έγγραφα/ Forsiktig! Se medfølgende dokumentasjon/ Atentie! Se va consulta documentația anexată/ Pozor! Prečítajte si sprievodnú dokumentáciu/ Varning! Se bifogad dokumentation/ Figyelmeztetés! Olvassa el a kísérődokumentációt/ Awasi! Sila rukuj kepada dokumen yang disertakan dengan peranti ini/ Анхааруула! Хабарсан бичиг баримтай танилц/ ข้อควรระวัง! โปรดศึกษาเอกสารที่แนบมาด้วย/ Thận trọng! Vui lòng tham khảo tài liệu kèm theo/ 注意! 請參考隨附文件/ 주의 함께 제공되는 설명서 참고

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
ENGLISH - Instructions for use

Read all instructions prior to use. To avoid complications, observe all warnings and precautions throughout these instructions. Caution: US Federal Law restricts this device to sale by or on the order of a physician. In accordance with European Regulation (EU) 2017/745 on medical devices, In the event of any serious incident occurring in relation to this medical device, report the details to the manufacturer and the competent authority of the European Member State in which the user is established.

Product Description

The SEQUIRE® microcatheters are single use microcatheters primarily comprised of a luer lock hub, a strain relief cover and tube, central shaft, and distal tip with a radiopaque marker for visualization. The SEQUIRE® microcatheters are used for infusion of therapeutic agents, embolic material and fluid such as contrast media.

The catheters consist of a metal braid reinforced shaft with external multi-polymer tubing. The braid provides trackability, pushability and torque transmission to the microcatheter as well as visibility under fluoroscopy. The inner lumen is made of PTFE (polytetrafluoroethylene), which allows for the smooth passage of fluids, embolic agents and devices such as guide wires. The distal section of the shaft is coated in a hydrophilic polymer layer, which ensures high lubricity when wet with saline or blood.

At the proximal end there is a luer-compatible hub connected to a strain relief cover which helps prevent kinking due to the weight of the hub. The SEQUIRE®'s distal end has two radiopaque markers that contribute to the fluoroscopic visualization of the distal tip of the microcatheter. The two radiopaque bands mark the boundaries of the SEQUIRE® side holes/slits.

The SEQUIRE® microcatheters are single lumen devices and are available in several different diameters and lengths..

Table 2: SEQUIRE® Microcatheter Specifications

Attribute	Specification								
Microcatheter OD (Proximal/Distal)	2.9/2.4 Fr (0.97/0.8 mm)			2.9/2.7 Fr (0.97/0.9 mm)			3.0/2.8 Fr (1.00/0.93 mm)		
Catheter Inner Diameter	0.022" (0.56 mm)			0.025" (0.64 mm)			0.027" (0.69 mm)		
Functional Length	105 cm	130 cm	150 cm	105 cm	130 cm	150 cm	105 cm	130 cm	150 cm
Maximal Compatible Guidewire OD	0.018" (0.46 mm)			0.021" (0.53 mm)			0.021" (0.53 mm)		
Recommended Guiding Catheter	Min. 0.038" (0.97 mm) guidewire compatible								
Compatible Embolic Microspheres OD	70-500 µm			70-500 µm			70-700 µm		
Compatible Embolic Coil Size	0.018"			0.018"			0.018"		

Indications For Use

The SEQUIRE® microcatheters are intended for the infusion of contrast media into all peripheral vessels. The SEQUIRE® microcatheters are also intended for drug infusion in intra-arterial therapy and infusion of embolic materials. The SEQUIRE® microcatheters should not be used in cerebral vessels.

Contraindications

Generally, angiography or intravascular therapy is contraindicated for, but not limited to, the patients listed below:

- Patients in the acute phase of myocardial infarction
- Patients with serious arrhythmia
- Patients with serious serum electrolyte imbalance
- Patients who in prior procedures have developed an adverse reaction to the injection of contrast media
- Patients with renal dysfunction
- Patients with coagulopathy or those whose blood has suffered a serious change in coagulation capability for some reasons
- Patients who cannot lie on their back on the operating table because of congestive heart failure or some respiratory disorder
- Patients with mental disease or those who are not expected to lie quietly during angiography
- Patients who are or may be pregnant
- Any other patients who are judged unsuitable for the procedure by the physician.
- Do not use in patients for whom anticoagulant and antiplatelet therapy is contraindicated.

Complications

Angiography or intravascular therapy may be accompanied by, but not limited to, the following:

- Headache
- Nausea and vomiting
- Fever and chill
- Abnormality in blood sampling tests
- Blood pressure drifting
- Shock
- Myocardial infarction
- Renal failure
- Infection and pain at the puncture site
- Inflammation with embolic material
- Cerebral edema
- Bradycardia
- Cerebral infarction from peripheral artery occlusion
- Hemorrhage, hematoma, arterio-venous fistula and/or false aneurysm at puncture site
- Spasm, artery perforation, aortic dissection and/or false aneurysm with the use of a guide wire or catheter
- Behavior disorder
- Death

Warnings

- Always flush the lumen of the guiding catheter and the microcatheter continuously with sterile heparinized saline solution prior to insertion. Residual contrast media or blood clots on the microcatheter surface reduce its lubricity, preventing smooth microcatheter movement. If flushing fails to restore surface lubricity, suspend the use of the microcatheter and remove it slowly and carefully, together with the guiding catheter.
- Do not pull the microcatheter with excessive force. Pulling the microcatheter with excessive force may cause serious mechanical damage, which may necessitate retrieval.
- Do not use excessive force when navigating the microcatheter or when advancing the guide wire through a kinked or blocked microcatheter. This may result in mechanical damage to the microcatheter and may cause damage to the blood vessels.
- Do not forcefully advance the microcatheter in extremely tortuous vasculature. This may result in kinks in the microcatheter or damage to the blood vessels.
- Manipulating the microcatheter and/or the guide wire against resistance may result in damage to the blood vessel, the microcatheter or the guide

wire. If the situation is not able to be resolved, withdraw the entire system of the microcatheter with the guiding catheter. During the procedure, continuously monitor the manipulation and location of the microcatheter in the vessel by confirming the position of the microcatheter tip through a high-resolution fluoroscope and/or a digital subtraction angiography monitor. If any resistance is felt, do not advance or withdraw the microcatheter until the cause of resistance is determined through a high-resolution fluoroscope and/or a digital subtraction angiography monitor.

- Do not soak or wipe the microcatheter with agents containing organic solvents, like antiseptic alcohol. It may damage the microcatheter or it may decrease the lubricity of the microcatheter.
- The infusion pressure in this microcatheter must never exceed 8270 kPa (1200 psi). A pressure that is higher than this maximum may cause serious mechanical failure in the microcatheter.
- The appropriate antiplatelet/anticoagulation therapy should be administered pre- and post-procedure in accordance with standard medical practice.

Precautions

- US Federal law restricts this device to sale by or on the order of a physician.
- The microcatheter must only be used by physicians who are thoroughly trained in the intended procedures.
- The entire procedure must be carried out in a sterile environment.
- Before opening, verify that the sterile packaging has not been compromised or damaged in any way.
- Before opening, ensure that the use-by date, printed on the label, has not expired.
- Before starting the procedure, the microcatheter itself must be visually inspected for any kinks, bends or other mechanical irregularities.
- Before use, the surface of the microcatheter must be completely wet with sterile heparinized saline solution to activate the hydrophilic coating.
- Manipulation of the microcatheter must be constantly monitored through a high-resolution fluoroscope and/or a digital subtraction angiography monitor.

- Refer to instructions for use information on any drugs and/or devices you may wish to use with this microcatheter to determine compatibility and to prevent microcatheter damage.
- Before starting a procedure, ensure that all devices and tools are in acceptable working condition.
- Choose the microcatheter size that is most suitable from a diagnostic and anatomical point of view.
- This device is sterilized with ethylene oxide gas (ETO) and is intended for single use only. Do not re-sterilize and/or reuse this device, as this might lead to patient infection or injury.
- Do not use if the unit package or the product has been damaged or soiled.
- Use immediately after opening the package and dispose of safely following your local procedure for the disposable of medical waste.

Precautions for storage

- Avoid exposure to water, direct sunlight, extreme temperature or humidity during storage.

Directions for use

1. Carefully remove the microcatheter in its hoop holder from the sterile pouch.
2. Flush the holder with sterile heparinized saline solution through the attached hub on the holder using a syringe, or immerse the microcatheter in its holder in a sterile heparinized saline solution bath to thoroughly wet the surface of the microcatheter.

CAUTION: The heparinized saline solution should be injected slowly so that the microcatheter is not driven out of the holder.

3. Slowly remove the microcatheter from its holder.

If any resistance is felt, do not use force. Inject sterile heparinized saline solution into the holder again and try to remove gently once more.

CAUTION: Do not use if the microcatheter has been damaged or if any other abnormality is observed. Hold the microcatheter by its hub during handling.

4. Use a syringe to prime the microcatheter lumen with sterile heparinized saline solution through its hub. To reduce the injection resistance, the use of a 1 ml or 2.5 ml luer lock syringes is recommended. Slowly inject 2-3 ml into the microcatheter until more than 10 drops of the solution appear out of its tip. Priming is completed if no air bubbles can be seen in these drops.
5. Attach a hemostatic valve or rotating hemostatic valve (Tuohy-Borst type) to the hub of the microcatheter, if necessary. Insert a guide wire, previously immersed in sterile heparinized saline solution and of a compatible size into the microcatheter through its hub or the attached valve and advance the guide wire to the distal end of the microcatheter. A torque device may be attached to the proximal end of the guide wire to facilitate guide wire manipulation. To maintain surface lubricity, immerse the microcatheter and guide wire assembly in a sterile heparinized saline solution bath or put it back into the microcatheter holder filled with sterile heparinized saline solution.

CAUTION: Do not insert the guide wire through the microcatheter's distal end. This may damage the microcatheter. When applicable, prime the hemostatic valve first, before inserting the guide wire into the microcatheter and advancing it to the distal end of the microcatheter.

6. Insert a guiding catheter into the patient's blood vessel. Attach a rotating hemostatic valve (Tuohy-Borst type) to the guiding catheter and irrigate the catheter continuously with heparinized saline solution. Antiplatelet/anticoagulation therapy should be administered in accordance with standard medical practice.

Insert the microcatheter and the guide wire assembly through the valve into the guiding catheter and advance to the distal end of the guiding catheter. For smooth insertion through the rotating hemostatic valve and guiding catheter, it is recommended to keep the guide wire tip within the microcatheter until the microcatheter reaches the distal end of the guiding catheter.

WARNING: Do not manipulate and/or withdraw the microcatheter through a metal entry needle or a metal dilator. Manipulation and/or withdrawal through a metal entry needle or a metal dilator may result in abrasion of the surface coating, destruction and/or separation of the microcatheter shaft. If the guiding catheter is fitted with a stopcock, do not close the stopcock with the microcatheter inside the guiding catheter. The microcatheter may be broken. Make sure that the guiding catheter does not slip out of the blood vessel. If the guiding catheter is dislocated from the vessel while the microcatheter and/or the guide wire are being advanced, this may result in damage to the microcatheter.

CAUTION: Do not excessively tighten the rotating hemostatic valve on the microcatheter, or manipulate the microcatheter through a tightened valve. Damage to the microcatheter may occur. If resistance is felt, do not force the microcatheter into the guiding catheter as this may result in the damage of

Directions for use (cont.)

the microcatheter. Do not advance the microcatheter without a guide wire, as kinks in the distal and proximal sections of the microcatheter may occur.

7. Continuously monitor the location of the microcatheter tip in the blood vessel during the procedure using high-resolution fluoroscope and/or a digital subtraction angiography monitor.

WARNING: If any resistance is felt, do not advance or withdraw the microcatheter until the cause of resistance is determined through a high-resolution fluoroscope and/or a digital subtraction angiography monitor. Manipulating the microcatheter and/or the guide wire against resistance may result in damage to the blood vessel, microcatheter or guide wire. If the microcatheter is advanced inside the blood vessel without the guide wire, it may result in damage to the vessel. When re-inserting the guide wire into the microcatheter, carefully advance the guide wire while monitoring the position of the guide wire tip through a high-resolution fluoroscope and/or a digital subtraction angiography monitor. Any erratic movement may cause damage to the vessel.

CAUTION: When advancing the microcatheter into the peripheral vessel, draw it back slightly under fluoroscopy each time it has been advanced, to ensure that the microcatheter has not been advanced so far that it cannot be withdrawn. Do not manipulate the microcatheter by force. The microcatheter tip is highly flexible and may be stretched or damaged.

8. When the desired site is reached, remove the guide wire from the microcatheter.

CAUTION: Drawing back the guide wire against resistance may cause the microcatheter to kink. If any resistance is felt, draw the microcatheter back to a position where no guide wire resistance is felt, then remove the guide wire. If the guide wire is withdrawn without this manipulation, the microcatheter may be impaired. Rinse residual blood from the removed guide wire in a sterile heparinized saline solution bath. If the residual stains will not come off, wipe the guide wire once with a sterile gauze moistened with heparinized saline solution. Blood remaining on the guide wire may cause resistance when re-inserted back into the microcatheter.

9. Before introducing an embolic material or other agent, slowly inject a small volume of contrast media into the microcatheter using a syringe and verify under high resolution fluoroscopy and/or a digital subtraction angiography monitor that the media is dispelled from the distal microcatheter tip. With its small lumen, the microcatheter offers a high resistance to infusion. When injecting contrast media or drugs with a syringe, use a 1 ml syringe or smaller. Refer to instructions for use for any drugs and/or devices you may wish to use with this microcatheter to determine compatibility and prevent microcatheter damage. If different types of embolic material or agents are being used, it is recommended to use a new microcatheter each time. When injecting contrast media through a SEQUIRE® microcatheter under fluoroscopy, it may be seen exiting radially from the slits which are located between the two radiopaque markers, in addition to exiting through the distal tip.

Note: It is recommended to heat the contrast medium to 37 °C before use.

WARNING: If an increase of resistance is felt when injecting fluid through the microcatheter, replace the microcatheter with a new one. Injection against increased resistance may cause the microcatheter to rupture, potentially

injuring the blood vessel. If no contrast media is visible flowing from the microcatheter, it may indicate kinking of the microcatheter. If drawing back the microcatheter fails to correct the kink, replace the microcatheter with a new one. Do not try to correct the kink by inserting a guide wire or by power injection. Introduction of embolic materials or agents without first correcting the kink, or attempts to correct the kink by inserting a guide wire or by power injection may cause the microcatheter to rupture or separate and this may result in damage to the vessel. Friction between the microcatheter's inner lumen and the embolic material may cause the microcatheter to advance forward, potentially resulting in perforation of the vessel wall. To prevent this, draw the microcatheter back slightly and hold in place.

CAUTION: Increased resistance to infusion suggests that the microcatheter may be blocked by the drug or contrast media being infused or with blood clots. Immediately terminate the infusion and replace the microcatheter. When a power injector is to be used, follow the instructions given below under "Instruction for Using a Power Injector with the microcatheter". In case of using organic solvents, make sure to check its compatibility by checking its labeling before use. Before use, ensure the compatibility of the size of the embolic microspheres to the chosen microcatheter size. The compatible microsphere size range is indicated on the microcatheter's external packaging label and in Table 2. Before use, confirm the size compatibility of the coiled embolic material and supportive device and the microcatheter. When introducing embolic coils or material, do not use devices or materials exceeding the maximal dimensions detailed in Table 2. Always check the movement of the embolic material and supportive device through a high-resolution fluoroscope and/or a digital subtraction angiography monitor. Do not advance or withdraw the microcatheter if any resistance is felt in the blood vessel, especially while using embolic material. Advance or withdraw the microcatheter, only after the cause of resistance is determined through a high-resolution fluoroscope and/or a digital subtraction angiography monitor. Any quick and unreasonable movement may cause the microcatheter to break/rupture/separate, which may result in damage to the vessel.

10. Before inserting the microcatheter into additional vessels, flush the microcatheter with sterile heparinized saline solution. If any resistance is felt during the insertion of the guide wire, do not advance the guide wire any farther and replace the microcatheter. If there is any difficulty in inserting the guide wire into the microcatheter hub, the guide wire and/or the hub can be rotated slightly to facilitate insertion.

WARNING: When re-inserting the guide wire into the microcatheter, verify the location of the guide wire tip through a high-resolution fluoroscope and/or digital subtraction angiography monitor. Any erratic movement of the wire may cause mechanical damage to the microcatheter, which may result in damage to the vessel.

11. When the procedure is completed, carefully remove the microcatheter together with the guiding catheter.

WARNING: If any resistance is felt, do not remove the microcatheter by force. Withdraw the microcatheter carefully together with the guiding catheter. Removing the microcatheter by force may result in the fracture of the microcatheter, which may necessitate retrieval.

Instructions for Using a Power Injector with the Microcatheter

A power injector can be used to infuse a contrast media through the microcatheter. Observe the warnings and cautions given below. The flow rate depends upon factors such as the viscosity of the contrast media, the type and temperature of the media, the model and setting of the power injector and how the injector is connected to the microcatheter.

WARNING: Do not use a power injector to infuse agents other than contrast media, as the microcatheter may become blocked. The injection pressure must not exceed the maximum injection pressure that corresponds to outer diameter of each microcatheter tip (see Table 3 below). Exceeding the maximal injection pressure may cause the microcatheter to rupture.

Table 3: Operational Information

Catheter OD (distal)	Usable Length of catheter (cm)	Maximum Injection Pressure	Dead Space Volume [ml]	Flow Rate with Iodinated Contrast Media (Iohexol 300 mg/ml, 37°C) [ml/sec]
2.4 Fr	105	8270 kPa (1200 psi)	0.48	2.7
	130		0.54	2.2
	150		0.59	1.9
2.7 Fr	105	8270 kPa (1200 psi)	0.56	3.8
	130		0.64	3.3
	150		0.7	3.1
2.8 Fr	105	8270 kPa (1200 psi)	0.62	4.8
	130		0.71	3.8
	150		0.78	3.5

- Under high-resolution fluoroscopy and DSA monitor, inject a small amount of contrast media with a syringe and confirm the flow of contrast media out of the microcatheter tip before using the power injector.
- If an expansion of the microcatheter outer diameter is observed during the injection, the pressure may have exceeded the maximum pressure limit. In such a case, stop the injection immediately.
- When securing the microcatheter in position, secure it by the hub so that microcatheter shaft is not damaged. When securing, do not hold the microcatheter shaft with forceps, as this may result in microcatheter separation.

CAUTION: If the microcatheter has been kinked or bent, it must be replaced.

- Connect the power injector to the microcatheter using a pressure-resistant extension tube.
- When re-inserting the guide wire after completion of angiography, the microcatheter’s inner lumen should be flushed with sterile heparinized saline solution.

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ACCURATE MEDICAL THERAPEUTICS LTD.
Distributed in Australia by Guerbet Australia Pty Ltd, Sydney NSW, Australia